

For Research Use Only

# CEP350-Specific Recombinant antibody, PBS Only (Capture)



Catalog Number: 83208-3-PBS

## Basic Information

Catalog Number:

83208-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_014810

GeneID (NCBI):

9857

UNIPROT ID:

Q5VT06

Full Name:

centrosomal protein 350kDa

Calculated MW:

351 kDa

Purification Method:

Protein A purification

CloneNo.:

230320E1

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00178-1, CEP350-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83208-3-PBS. Detection antibody: 83208-2-PBS. Standard: SY00642. Range: 0.313-20 ng/mL.