

For Research Use Only

Mouse Il6ra Recombinant antibody, PBS Only (Detector)

Catalog Number: 83209-1-PBS



Basic Information

Catalog Number:

83209-1-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_010559.3

GeneID (NCBI):

16194

UNIPROT ID:

P22272

Full Name:

interleukin 6 receptor, alpha

Calculated MW:

50KD

Purification Method:

Protein A purification

CloneNo.:

23037188

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00182-1, Mouse CD126/IL-6R alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83209-2-PBS. Detection antibody: 83209-1-PBS. Standard: Eg0664. Range: 1.25-40 ng/mL