

For Research Use Only

IGDCC3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83211-1-PBS

Basic Information

Catalog Number: 83211-1-PBS	GenBank Accession Number: BC067107	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9543	CloneNo.: 23044988
Source: Rabbit	UNIPROT ID: Q8IVU1	
Isotype: IgG	Full Name: immunoglobulin superfamily, DCC subclass, member 3	
Immunogen Catalog Number: AG10276	Calculated MW: 814 aa, 87 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00192-1, IGDC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83211-2-PBS. Detection antibody: 83211-1-PBS. Standard: Ag10276. Range: 1.25-20 ng/mL



Cytometric bead array standard curve ofMP00192-2, IGDC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83211-3-PBS. Detection antibody: 83211-1-PBS. Standard: Ag10276. Range: 1.25-20 ng/mL