

For Research Use Only

CBP80 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83212-1-PBS



Basic Information

Catalog Number:

83212-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0496

GenBank Accession Number:

BC001450

GeneID (NCBI):

4686

UNIPROT ID:

Q09161

Full Name:

nuclear cap binding protein subunit 1,

80kDa

Calculated MW:

92 kDa

Purification Method:

Protein A purification

CloneNo.:

240012A4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00251-1, CBP80 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83212-3-PBS. Detection antibody: 83212-1-PBS. Standard: Ag0496. Range: 2.5-20 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83212-1-PBS against Human CBP80 was performed. The affinity constant is below 1 pM.