

For Research Use Only

Mouse CD86 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83213-2-PBS

Basic Information

Catalog Number:

83213-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_019388

GeneID (NCBI):

12524

UNIPROT ID:

P42082-1

Full Name:

CD86 antigen

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

230476G6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00195-2, Mouse Cd86 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83213-2-PBS. Detection antibody: 83213-3-PBS. Standard: Eg0802. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00195-3, Mouse Cd86 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83213-2-PBS. Detection antibody: 83213-1-PBS. Standard: Eg0802. Range: 0.625-20 ng/mL