

For Research Use Only

Mouse CD86 Recombinant monoclonal antibody, PBS Only

Catalog Number: 83213-3-PBS



Basic Information

Catalog Number:

83213-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0802

GenBank Accession Number:

NM_019388

GeneID (NCBI):

12524

UNIPROT ID:

P42082-1

Full Name:

CD86 antigen

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

230476D8

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

