

For Research Use Only

Mouse B7-2/CD86 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 83213-9-PBS

Basic Information

Catalog Number:

83213-9-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0802

GenBank Accession Number:

NM_019388

GeneID (NCBI):

12524

UNIPROT ID:

P42082-1

Full Name:

CD86 antigen

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

251860F8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Serum of six mice was measured. The CD86 concentration of detected samples was determined to be 3,467.1 pg/mL with a range of 550.4-7,806.3 pg/mL.



Sandwich ELISA standard curve of MP00195-5, Mouse CD86 Recombinant Matched Antibody Pair - PBS only. 83213-10-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0802. 83213-9-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.