

For Research Use Only

# PDGF-BB Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 83214-3-PBS



## Basic Information

Catalog Number:

83214-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25986

GenBank Accession Number:

BC029822

GeneID (NCBI):

5155

UNIPROT ID:

P01127-1

Full Name:

platelet-derived growth factor beta polypeptide (simian sarcoma viral (v-sis) oncogene homolog)

Calculated MW:

241 aa, 27 kDa

Purification Method:

Protein A purification

CloneNo.:

230486C9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00198-1, PDGF-BB Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83214-2-PBS. Detection antibody: 83214-3-PBS. Standard: Ag25986. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00198-2, PDGF-BB Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83214-4-PBS. Detection antibody: 83214-3-PBS. Standard: Ag25986. Range: 0.625-20 ng/mL