

For Research Use Only

ZFYVE26 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83215-1-PBS

Basic Information

Catalog Number:

83215-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_015346

GeneID (NCBI):

23503

UNIPROT ID:

Q68DK2

Full Name:

zinc finger, FYVE domain containing 26

Calculated MW:

285 kDa

Purification Method:

Protein A purification

CloneNo.:

230493F2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00202-1, ZFYVE26 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83215-2-PBS. Detection antibody: 83215-1-PBS. Standard: SY01497. Range: 2.5-40 ng/mL