

For Research Use Only

# NDUFB8 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83216-2-PBS



## Basic Information

Catalog Number:

83216-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6495

GenBank Accession Number:

BC000466

GeneID (NCBI):

4714

UNIPROT ID:

O95169

Full Name:

NADH dehydrogenase (ubiquinone) 1  
beta subcomplex, 8, 19kDa

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

230510F6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00204-1, NDUF88 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83216-1-PBS. Detection antibody: 83216-2-PBS. Standard: Ag6495. Range: 2.5-20 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83216-2-PBS against Human NDUF88 was performed. The affinity constant is below 1 pM.