

For Research Use Only

DNAJC11 Recombinant antibody, PBS Only (Detector)

Catalog Number:83218-5-PBS



Basic Information

Catalog Number:

83218-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11310

GenBank Accession Number:

BC014145

GeneID (NCBI):

55735

UNIPROT ID:

Q9NVH1

Full Name:

DnaJ (Hsp40) homolog, subfamily C, member 11

Calculated MW:

507 aa, 57 kDa

Purification Method:

Protein A purification

CloneNo.:

230492A1

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01315-3, Human DNAJC11 Recombinant Matched Antibody Pair - PBS only. 83218-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag11310. 83218-5-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL



HepG2 cell lysate was measured. The human DNAJC11 concentration of detected samples was determined to be 24.80 ng/mL (based on a 1.30 mg/mL extract load).