

For Research Use Only

nectin 4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83223-2-PBS



Basic Information

Catalog Number:

83223-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC010423

GeneID (NCBI):

81607

UNIPROT ID:

Q96NY8

Full Name:

poliovirus receptor-related 4

Calculated MW:

510 aa, 55 kDa

Purification Method:

Protein A purification

CloneNo.:

240009G5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Nectin-4, also known as PVRL4, is a type-I transmembrane glycoprotein that belongs to the nectin subfamily of immunoglobulin-like adhesion molecules that participate in Ca^{2+} -independent cell-cell adhesion (PMID: 18703497). The extracellular domain of Nectin-4, which contains three immunoglobulin-like domains (V and two C2-type domains, VCC), can be proteolytically cleaved to release a soluble form. Nectin-4 interacts with afadin via its carboxyl-terminal cytoplasmic sequence and trans-interacts with nectin-1/PRR1 through V domain interaction (PMID: 11544254). It acts as a receptor for measles virus (PMID: 22048310). Nectin-4 is overexpressed in multiple human malignancies and the aberrant expression is correlated with cancer progression and poor prognosis (PMID: 35131876). Nectin-4 can be detected as a soluble form of 43 kDa and a transmembrane form of 66 kDa (PMID: 15784625).

Storage

Storage:

Store at -80°C .

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00147-3, Nectin-4/PVRL4 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83223-5-PBS. Detection antibody: 83223-2-PBS. Standard: Eg0415. Range: 1.25-20 ng/mL