

For Research Use Only

PARK2/Parkin Recombinant antibody, PBS Only (Detector)

Catalog Number: 83224-1-PBS



Basic Information

Catalog Number:

83224-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5179

GenBank Accession Number:

BC022014

GeneID (NCBI):

5071

UNIPROT ID:

O60260

Full Name:

Parkinson disease (autosomal
recessive, juvenile) 2, parkin

Calculated MW:

52 kDa

Purification Method:

Protein A purification

CloneNo.:

240077B8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00169-1, Parkin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83224-3-PBS. Detection antibody: 83224-1-PBS. Standard: Ag5179. Range: 5-80 ng/mL



Cytometric bead array standard curve of MP00169-2, Parkin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83224-2-PBS. Detection antibody: 83224-1-PBS. Standard: Ag5179. Range: 5-80 ng/mL