

For Research Use Only

# PARK2/Parkin Recombinant antibody, PBS Only (Capture)

Catalog Number: 83224-3-PBS



## Basic Information

Catalog Number:

83224-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5179

GenBank Accession Number:

BC022014

GeneID (NCBI):

5071

UNIPROT ID:

O60260

Full Name:

Parkinson disease (autosomal  
recessive, juvenile) 2, parkin

Calculated MW:

52 kDa

Purification Method:

Protein A purification

CloneNo.:

240077G2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00169-1, Parkin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83224-3-PBS. Detection antibody: 83224-1-PBS. Standard: Ag5179. Range: 5-80 ng/mL