

For Research Use Only

SULT1B1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83228-2-PBS



Basic Information

Catalog Number:

83228-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9009

GenBank Accession Number:

BC010895

GeneID (NCBI):

27284

UNIPROT ID:

O43704

Full Name:

sulfotransferase family, cytosolic, 1B, member 1

Calculated MW:

296 aa, 35 kDa

Purification Method:

Protein A purification

CloneNo.:

240119F6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00184-1, SULT1B1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83228-1-PBS. Detection antibody: 83228-2-PBS. Standard: Ag9009. Range: 2.5-80 ng/mL



Cytometric bead array standard curve of MP00184-2, SULT1B1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83228-2-PBS. Detection antibody: 83228-3-PBS. Standard: Ag9009. Range: 2.5-80 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83228-2-PBS against Human SULT1B1 was performed. The affinity constant is 5.92 nM.