

For Research Use Only

SP8 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83230-4-PBS



Basic Information

Catalog Number:

83230-4-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22934

GenBank Accession Number:

BC038669

GeneID (NCBI):

221833

UNIPROT ID:

Q8IXZ3

Full Name:

Sp8 transcription factor

Calculated MW:

508 aa, 51 kDa

Purification Method:

Protein A purification

CloneNo.:

240011B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00190-2, SP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83230-4-PBS. Detection antibody: 83230-1-PBS. Standard: Ag22934. Range: 2.5-20 ng/mL.