

For Research Use Only

# GBX1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83232-2-PBS



## Basic Information

Catalog Number:

83232-2-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001098834

GeneID (NCBI):

2636

UNIPROT ID:

Q14549

Full Name:

gastrulation brain homeobox 1

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

240108H1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00201-1, GBX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83232-5-PBS. Detection antibody: 83232-2-PBS. Standard: SY01055. Range: 0.313-5 ng/mL.



Cytometric bead array standard curve of MP00201-3, GBX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83232-2-PBS. Detection antibody: 83232-1-PBS. Standard: SY01055. Range: 1.25-20 ng/mL.