

For Research Use Only

# RELN Recombinant antibody, PBS Only (Detector)

Catalog Number: 83233-2-PBS



## Basic Information

Catalog Number:

83233-2-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_005045

GeneID (NCBI):

5649

UNIPROT ID:

P78509

Full Name:

reelin

Calculated MW:

388 kDa

Purification Method:

Protein A purification

CloneNo.:

230387D10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00207-1, RELN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83233-1-PBS. Detection antibody: 83233-2-PBS. Standard: SY00905. Range: 1.25-20 ng/mL.