

For Research Use Only

# MKNK2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83235-3-PBS



## Basic Information

Catalog Number:

83235-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10546

GenBank Accession Number:

BC018345

GeneID (NCBI):

2872

UNIPROT ID:

Q9HBH9

Full Name:

MAP kinase interacting  
serine/threonine kinase 2

Calculated MW:

465 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

240045B9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00216-1, MKNK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83235-4-PBS. Detection antibody: 83235-3-PBS. Standard: Ag10546. Range: 1.25-40 ng/mL.



Cytometric bead array standard curve of MP00216-3, MKNK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83235-3-PBS. Detection antibody: 83235-2-PBS. Standard: Ag10546. Range: 1.25-40 ng/mL.