

For Research Use Only

TDO2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83236-3-PBS

Basic Information

Catalog Number: 83236-3-PBS	GenBank Accession Number: BC005355	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 6999	CloneNo.: 240068G7
Source: Rabbit	UNIPROT ID: P48775	
Isotype: IgG	Full Name: tryptophan 2,3-dioxygenase	
Immunogen Catalog Number: AG8531	Calculated MW: 406 aa, 48 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00219-2, TDO2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83236-2-PBS. Detection antibody: 83236-3-PBS. Standard: Ag8531. Range: 1.25-40 ng/mL.