

For Research Use Only

SLC22A10 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83237-1-PBS

Basic Information

Catalog Number: 83237-1-PBS	GenBank Accession Number: NM_001039752	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 387775	CloneNo.: 240083A3
Source: Rabbit	UNIPROT ID: Q63ZE4	
Isotype: IgG	Full Name: solute carrier family 22, member 10	
Immunogen Catalog Number: AG34521	Calculated MW: 541 aa, 60 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00221-1, SLC22A10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83237-3-PBS. Detection antibody: 83237-1-PBS. Standard: Ag34521. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00221-2, SLC22A10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83237-1-PBS. Detection antibody: 83237-2-PBS. Standard: Ag34521. Range: 0.625-20 ng/mL