

For Research Use Only

# SLC22A10 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83237-3-PBS

## Basic Information

Catalog Number:

83237-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34521

GenBank Accession Number:

NM\_001039752

GeneID (NCBI):

387775

UNIPROT ID:

Q63ZE4

Full Name:

solute carrier family 22, member 10

Calculated MW:

541 aa, 60 kDa

Purification Method:

Protein A purification

CloneNo.:

240083F4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve ofMP00221-1, SLC22A10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83237-3-PBS. Detection antibody: 83237-1-PBS. Standard: Ag34521. Range: 0.625-20 ng/mL.