

For Research Use Only

SOX8 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83238-1-PBS

Basic Information

| | | |
|--------------------------------------|--|--|
| Catalog Number: 83238-1-PBS | GenBank Accession Number: BC031797 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 30812 | CloneNo.: 240120G2 |
| Source: Rabbit | UNIPROT ID: P57073 | |
| Isotype: IgG | Full Name: SRY (sex determining region Y)-box 8 | |
| Immunogen Catalog Number: AG34395 | Calculated MW: 446 aa, 47 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00222-1, SOX8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83238-1-PBS. Detection antibody: 83238-3-PBS. Standard: Ag34395. Range: 0.625-20 ng/mL