

For Research Use Only

SOX8 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83238-2-PBS

Basic Information

Catalog Number: 83238-2-PBS	GenBank Accession Number: BC031797	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 30812	CloneNo.: 240120H5
Source: Rabbit	UNIPROT ID: P57073	
Isotype: IgG	Full Name: SRY (sex determining region Y)-box 8	
Immunogen Catalog Number: AG34395	Calculated MW: 446 aa, 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00222-2, SOX8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83238-3-PBS. Detection antibody: 83238-2-PBS. Standard: Ag34395. Range: 0.625-20 ng/mL



Biolayer interferometry (BLI) kinetic assays of 83238-2-RR against Human SOX8 were performed. The affinity constant is 6.54 nM.