

For Research Use Only

# SOX8 Recombinant antibody

Catalog Number: 83238-2-RR



## Basic Information

Catalog Number:

83238-2-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34395

GenBank Accession Number:

BC031797

GeneID (NCBI):

30812

UNIPROT ID:

P57073

Full Name:

SRX (sex determining region Y)-box 8

Calculated MW:

446 aa, 47 kDa

Purification Method:

Protein A purification

CloneNo.:

240120H5

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83238-2-RR against Human SOX8 were performed. The affinity constant is 6.54 nM.