

For Research Use Only

# SLC16A11 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83239-3-PBS

## Basic Information

Catalog Number: 83239-3-PBS	GenBank Accession Number: BC093860	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 162515	CloneNo.: 240081C8
Source: Rabbit	UNIPROT ID: Q8NCK7	
Isotype: IgG	Full Name: solute carrier family 16, member 11 (monocarboxylic acid transporter 11)	
Immunogen Catalog Number: AG24920	Calculated MW: 471 aa, 48 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve ofMP00223-1, SLC16A11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83239-3-PBS. Detection antibody: 83239-2-PBS. Standard: Ag24920. Range: 0.625-20 ng/mL.