

For Research Use Only

# SLC16A11 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83239-4-PBS

## Basic Information

Catalog Number:

83239-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24920

GenBank Accession Number:

BC093860

GeneID (NCBI):

162515

UNIPROT ID:

Q8NCK7

Full Name:

solute carrier family 16, member 11  
(monocarboxylic acid transporter 11)

Calculated MW:

471 aa, 48 kDa

Purification Method:

Protein A purification

CloneNo.:

240081H8

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve ofMP00223-3, SLC16A11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83239-4-PBS. Detection antibody: 83239-2-PBS. Standard: Ag24920. Range: 1.25-20 ng/mL