

For Research Use Only

SLC16A11 Recombinant antibody, PBS Only



Catalog Number: 83239-5-PBS

Basic Information

Catalog Number:

83239-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24920

GenBank Accession Number:

BC093860

GeneID (NCBI):

162515

UNIPROT ID:

Q8NCK7

Full Name:

solute carrier family 16, member 11
(monocarboxylic acid transporter 11)

Calculated MW:

471 aa, 48 kDa

Purification Method:

Protein A purification

CloneNo.:

240081H5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83239-5-RR (SLC16A11 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83239-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83239-5-RR against Human SLC16A11 were performed. The affinity constant is below 1 pM.