

For Research Use Only

Glypican 2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83241-2-PBS



Basic Information

Catalog Number:

83241-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4041

GenBank Accession Number:

BC027972

GeneID (NCBI):

221914

UNIPROT ID:

Q8N158

Full Name:

glypican 2

Calculated MW:

579 aa, 63 kDa

Purification Method:

Protein A purification

CloneNo.:

230518B1

Affinity:

$K_D < 1.00 \times 10^{-12}$

$K_{off} < 1.00 \times 10^{-6}$

$K_{on} = 4.45 \times 10^4$

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00156-2, Glypican 2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83241-1-PBS. Detection antibody: 83241-2-PBS. Standard: Ag4041. Range: 2.5-80 ng/mL