

For Research Use Only

MLLT1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83242-1-PBS

Basic Information

Catalog Number:

83242-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34210

GenBank Accession Number:

BC148806

GeneID (NCBI):

4298

UNIPROT ID:

Q03111

Full Name:

myeloid/lymphoid or mixed-lineage

leukemia (trithorax homolog,

Drosophila); translocated to, 1

Calculated MW:

62 kDa

Purification Method:

Protein A purification

CloneNo.:

230519B12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00158-3, MLLT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83242-1-PBS. Detection antibody: 83242-2-PBS. Standard: Ag34210. Range: 1.25-40 ng/mL.