

For Research Use Only

MLLT1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83242-2-PBS

Basic Information

Catalog Number: 83242-2-PBS	GenBank Accession Number: BC148806	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 4298	CloneNo.: 230519F3
Source: Rabbit	UNIPROT ID: Q03111	
Isotype: IgG	Full Name: myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila); translocated to, 1	
Immunogen Catalog Number: AG34210	Calculated MW: 62 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00158-1, MLLT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83242-3-PBS. Detection antibody: 83242-2-PBS. Standard: Ag34210. Range: 1.25-40 ng/mL



Cytometric bead array standard curve ofMP00158-3, MLLT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83242-1-PBS. Detection antibody: 83242-2-PBS. Standard: Ag34210. Range: 1.25-40 ng/mL