

For Research Use Only

Mouse Cxcl15 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83245-2-PBS



Basic Information

Catalog Number:

83245-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33964

GenBank Accession Number:

BC061138

GeneID (NCBI):

20309

UNIPROT ID:

Q9WVL7

Full Name:

chemokine (C-X-C motif) ligand 15

Calculated MW:

19 kDa

Purification Method:

Protein A purification

CloneNo.:

240163F1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00193-1, Mouse Cxcl15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83245-3-PBS. Detection antibody: 83245-2-PBS. Standard: Ag33964. Range: 0.625-20 ng/mL.