

For Research Use Only

# Mouse Cxcl15 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83245-3-PBS



## Basic Information

Catalog Number:

83245-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33964

GenBank Accession Number:

BC061138

GeneID (NCBI):

20309

UNIPROT ID:

Q9WVL7

Full Name:

chemokine (C-X-C motif) ligand 15

Calculated MW:

19 kDa

Purification Method:

Protein A purification

CloneNo.:

240163H5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00193-1, Mouse Cxcl15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83245-3-PBS. Detection antibody: 83245-2-PBS. Standard: Ag33964. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00193-2, Mouse Cxcl15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83245-3-PBS. Detection antibody: 83245-1-PBS. Standard: Ag33964. Range: 1.25-20 ng/mL