

For Research Use Only

# APCDD1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83246-3-PBS



## Basic Information

Catalog Number:

83246-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10303

GenBank Accession Number:

BC053324

GeneID (NCBI):

147495

UNIPROT ID:

Q8J025

Full Name:

adenomatosis polyposis coli down-regulated 1

Calculated MW:

514 aa, 59 kDa

Purification Method:

Protein A purification

CloneNo.:

240085C8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00199-1, APCDD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83246-3-PBS. Detection antibody: 83246-2-PBS. Standard: Ag10303. Range: 3.13-100 ng/mL.