

For Research Use Only

GJA5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83247-2-PBS

Basic Information

Catalog Number:

83247-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24199

GenBank Accession Number:

BC013313

GeneID (NCBI):

2702

UNIPROT ID:

P36382

Full Name:

gap junction protein, alpha 5, 40kDa

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

240072E10

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00208-1, GJA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83247-3-PBS. Detection antibody: 83247-2-PBS. Standard: Ag24199. Range: 6.25-100 ng/mL