

For Research Use Only

ABHD6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83249-1-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
83249-1-PBS	BC001698	Protein A purification
Size:	GeneID (NCBI):	CloneNo.:
1mg/ml	57406	230472C1
Source:	UNIPROT ID:	
Rabbit	Q9BV23	
Isotype:	Full Name:	
IgG	abhydrolase domain containing 6	
Immunogen Catalog Number:	Calculated MW:	
AG13968	337 aa, 38 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00225-1, ABHD6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83249-1-PBS. Detection antibody: 83249-2-PBS. Standard: Ag13968. Range: 3.125-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83249-1-PBS against Human ABHD6 was performed. The affinity constant is 0.384 nM.