

For Research Use Only

Ribophorin I Recombinant antibody, PBS Only (Detector)

Catalog Number: 83250-3-PBS



Basic Information

Catalog Number:

83250-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3576

GenBank Accession Number:

BC010839

GeneID (NCBI):

6184

UNIPROT ID:

P04843

Full Name:

ribophorin I

Calculated MW:

69 kDa

Purification Method:

Protein A purification

CloneNo.:

240044F12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00231-1, RPN1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83250-1-PBS. Detection antibody: 83250-3-PBS. Standard: Ag3576. Range: 2.5-80 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83250-3-PBS against Human Ribophorin I was performed. The affinity constant is 1.56 nM.