

For Research Use Only

NRG1, isoform Alpha Recombinant antibody, PBS Only (Detector)

Catalog Number: 83251-1-PBS



Basic Information

Catalog Number:

83251-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35114

GenBank Accession Number:

BC007675

GeneID (NCBI):

3084

UNIPROT ID:

Q02297

Full Name:

neuregulin 1

Calculated MW:

70 kDa

Purification Method:

Protein A purification

CloneNo.:

240067A4

Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP00232-6, Human NRG1, isoform Alpha Recombinant Matched Antibody Pair - PBS only. 83251-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag35114. 83251-1-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL



The mean NRG1, isoform Alpha concentration was determined to be 2,318.1 pg/mL in A549 cell extract based on a 2.5 mg/mL extract load.



Cytometric bead array standard curve of MP00232-2, NRG1, isoform Alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83251-3-PBS. Detection antibody: 83251-1-PBS. Standard: Ag35114. Range: 1.25-20 ng/mL