

For Research Use Only

ND6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83252-1-PBS

Basic Information

Catalog Number: 83252-1-PBS	GenBank Accession Number: YP_003024037	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 4541	CloneNo.: 240146A3
Source: Rabbit	UNIPROT ID: P03923	
Isotype: IgG	Full Name: NADH dehydrogenase, subunit 6 (complex I)	
Immunogen Catalog Number: AG34496	Calculated MW: 19 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00246-1, ND6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83252-1-PBS. Detection antibody: 83252-3-PBS. Standard: Ag34496. Range: 2.5-80 ng/mL.