

For Research Use Only

# ND6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83252-4-PBS

## Basic Information

Catalog Number:	83252-4-PBS	GenBank Accession Number:	YP_003024037	Purification Method:	Protein A purification
Size:	1mg/ml	GeneID (NCBI):	4541	CloneNo.:	240146F5
Source:	Rabbit	UNIPROT ID:	P03923		
Isotype:	IgG	Full Name:	NADH dehydrogenase, subunit 6 (complex I)		
Immunogen Catalog Number:	AG34496	Calculated MW:	19 kDa		

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00246-2, ND6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83252-4-PBS. Detection antibody: 83252-2-PBS. Standard: Ag34496. Range: 2.5-40 ng/mL.