

For Research Use Only

VSIG2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83255-1-PBS

Basic Information

Catalog Number: 83255-1-PBS	GenBank Accession Number: BC007313	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 23584	CloneNo.: 240139A2
Source: Rabbit	UNIPROT ID: Q96IQ7	Affinity: $K_D < 1.00 \times 10^{-12}$
Isotype: IgG	Full Name: V-set and immunoglobulin domain containing 2	$K_{off} < 1.00 \times 10^{-6}$ $K_{on} = 4.75 \times 10^4$
Immunogen Catalog Number: AG14342	Calculated MW: 327 aa, 34 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00161-1, VSG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83255-2-PBS. Detection antibody: 83255-1-PBS. Standard: Ag14342. Range: 1.25-20 ng/mL.