

For Research Use Only

TMEM40 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83256-3-PBS



Basic Information

Catalog Number:

83256-3-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27882

GenBank Accession Number:

BC020658

GeneID (NCBI):

55287

UNIPROT ID:

Q8WWA1

Full Name:

transmembrane protein 40

Purification Method:

Protein A purification

CloneNo.:

240102E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00163-1, TMEM40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83256-3-PBS. Detection antibody: 83256-1-PBS. Standard: Ag27882. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00163-2, TMEM40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83256-2-PBS. Detection antibody: 83256-3-PBS. Standard: Ag27882. Range: 0.625-20 ng/mL