

For Research Use Only

KCNJ15 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83257-2-PBS

Basic Information

| | | |
|--------------------------------------|--|--|
| Catalog Number: 83257-2-PBS | GenBank Accession Number: BC013327 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 3772 | CloneNo.: 240141C1 |
| Source: Rabbit | UNIPROT ID: Q99712 | |
| Isotype: IgG | Full Name: potassium inwardly-rectifying channel, subfamily J, member 15 | |
| Immunogen Catalog Number: AG34655 | Calculated MW: 375 aa, 43 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00164-2, KCNJ15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83257-1-PBS. Detection antibody: 83257-2-PBS. Standard: Ag34655. Range: 2.5-80 ng/mL.