

For Research Use Only

CBX5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83258-4-PBS



Basic Information

Catalog Number:

83258-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34810

GenBank Accession Number:

BC006821

GeneID (NCBI):

23468

UNIPROT ID:

P45973

Full Name:

chromobox homolog 5 (HP1 alpha

homolog, Drosophila)

Calculated MW:

191 aa, 22 kDa

Purification Method:

Protein A purification

CloneNo.:

240145D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00181-1, CBX5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83258-2-PBS. Detection antibody: 83258-4-PBS. Standard: Ag34810. Range: 3.125-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83258-4-PBS against Human CBX1 was performed. The affinity constant is below 1 pM.