

For Research Use Only

# SNAP25 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83259-3-PBS

## Basic Information

Catalog Number:	83259-3-PBS	GenBank Accession Number:	BC010647	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	6616	CloneNo.:	240069A3
Source:	Rabbit	UNIPROT ID:	P60880		
Isotype:	IgG	Full Name:	synaptosomal-associated protein,		
Immunogen Catalog Number:	AG6695	Calculated MW:	25kDa		
			23 kDa		

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00196-1, SNAP25 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83259-3-PBS. Detection antibody: 83259-1-PBS. Standard: Ag6695. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00196-3, SNAP25 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83259-3-PBS. Detection antibody: 83259-2-PBS. Standard: Ag6695. Range: 0.625-20 ng/mL