

For Research Use Only

PLK1-phospho Recombinant antibody, PBS Only (Capture)

Catalog Number: 83260-2-PBS



Basic Information

Catalog Number:

83260-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005030

GeneID (NCBI):

5347

UNIPROT ID:

P53350

Full Name:

polo-like kinase 1 (Drosophila)

Calculated MW:

68 kDa

Purification Method:

Protein A purification

CloneNo.:

240092D5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00238-1, Phospho-PLK1 (Ser326) Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83260-2-PBS. Detection antibody: 83260-3-PBS. Standard: SY00401. Range: 0.156-10 ng/mL



Cytometric bead array standard curve of MP00238-2, Phospho-PLK1 (Ser326) Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83260-2-PBS. Detection antibody: 83260-4-PBS. Standard: SY00401. Range: 0.313-10 ng/mL