

For Research Use Only

GNA12 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83263-2-PBS

Basic Information

Catalog Number: 83263-2-PBS	GenBank Accession Number: BC087537	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2768	CloneNo.: 240050D5
Source: Rabbit	UNIPROT ID: Q03113	
Isotype: IgG	Full Name: guanine nucleotide binding protein (G protein) alpha 12	
Immunogen Catalog Number: AG35068		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00276-2, GNA12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83263-1-PBS. Detection antibody: 83263-2-PBS. Standard: Ag35068. Range: 0.625-20 ng/mL