

For Research Use Only

DUS2L Recombinant antibody, PBS Only (Capture)

Catalog Number: 83265-1-PBS



Basic Information

Catalog Number:

83265-1-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23831

GenBank Accession Number:

BC006527

GeneID (NCBI):

54920

UNIPROT ID:

Q9NX74

Full Name:

dihydrouridine synthase 2-like, SMM1 homolog (S. cerevisiae)

Purification Method:

Protein A purification

CloneNo.:

240192A6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00261-2, DUS2L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83265-3-PBS. Detection antibody: 83265-4-PBS. Standard:Ag23831. Range: 1.25-40 ng/mL