

For Research Use Only

BUB3 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83266-3-PBS

Basic Information

Catalog Number:

83266-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25608

GenBank Accession Number:

BC005138

GeneID (NCBI):

9184

UNIPROT ID:

O43684

Full Name:

budding uninhibited by
benzimidazoles 3 homolog (yeast)

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

240087D6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00297-1, BUB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83266-3-PBS. Detection antibody: 83266-4-PBS. Standard: Ag25608. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00297-2, BUB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83266-1-PBS. Detection antibody: 83266-3-PBS. Standard: Ag25608. Range: 1.25-20 ng/mL