

For Research Use Only

# BUB3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83266-4-PBS

## Basic Information

Catalog Number: 83266-4-PBS	GenBank Accession Number: BC005138	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9184	CloneNo.: 240087F3
Source: Rabbit	UNIPROT ID: O43684	
Isotype: IgG	Full Name: budding uninhibited by benzimidazoles 3 homolog (yeast)	
Immunogen Catalog Number: AG25608	Calculated MW: 37 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve ofMP00297-1, BUB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83266-3-PBS. Detection antibody: 83266-4-PBS. Standard: Ag25608. Range: 0.625-20 ng/mL